

## Product Data Sheet

Product Name: DOP-DEDA

Cat. No.: GF13874

### Chemical Properties

Cas. No. 2247753-10-0

Formula C43H83N2O8P

M.Wt 787.12

Solubility

Storage Store at 2-8°C

**General tips** For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

**Shipping Condition** Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA