

Product Data Sheet

Product Name: Yaddle1
Cat. No.: GC74002

Chemical Properties

Cas. No.


Formula $C_{14}H_8ClF_3N_4OS$ M.Wt 372.75

Solubility DMSO : 100 mg/mL (268.28 mM; Need ultrasonic) Storage $-20^{\circ}C$

General tips For obtaining a higher solubility , please warm the tube at $37^{\circ}C$ and shake it in the ultrasonic bath for a while. Stock solution can be stored below $-20^{\circ}C$ for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Yaddle1 (compound 34) a novel Piezo1 agonist, and induces Ca^{2+} influx in human CD4^{+} T cell. Yaddle1 can be used for study of vaccine adjuvant.

References:

[1]. Goon S, et al. Exploring the Structural Attributes of Yoda1 for the Development of New-Generation Piezo1 Agonist Yaddle1 as a Vaccine Adjuvant Targeting Optimal T Cell Activation. J Med Chem. 2024;67(10):8225-8246.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA