

Product Data Sheet

Product Name: 1,1'-Bis(diphenylphosphino)ferrocene-palladium dichloride dichloromethane adduct
Cat. No.: GD02334

Chemical Properties

Cas. No. 95464-05-4

Formula

M.Wt


Solubility DMSO : 12.5 mg/mL (15.31 mM;
ultrasonic and warming and heat to
60°C)

Storage argon filling

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

Dichloro[1,1'-bis(diphenylphosphino)ferrocene]palladium(II) dichloromethane adduct is a biochemical reagent that can be used as a biological material or organic compound for life science related research.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA