

Product Data Sheet

Product Name: Perylene-3-boronic acid

Cat. No.: GF04796

Chemical Properties

Cas. No. 955121-20-7

Formula C₂₀H₁₃BO₂

M.Wt 296.12702


Solubility

Storage 2-8°C, away from moisture, sealed storage

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA