

Product Data Sheet

Product Name: KL044
Cat. No.: GC63033

Chemical Properties

Cas. No. 1801856-93-8


Formula $C_{21}H_{14}ClN_3O$ M.Wt 359.81

Solubility DMSO : 100 mg/mL (277.92 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

KL044, a stabilizer of the clock protein cryptochrome (CRY) , is a potent chemical probe with a pEC50 value of 7.32, leading to the extension of the circadian period and repression of Per2

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

activity[1].

KL044 (0-3.7 μ M; 8 h) effectively stabilizes the CRY1-LUC fusion protein without affecting the stability of LUC in the HEK293 stable cell line[1].

[1]. Jae Wook Lee, et al. Development of Small-Molecule Cryptochrome Stabilizer Derivatives as Modulators of the Circadian Clock. ChemMedChem. 2015 Sep;10(9):1489-97.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA