

Product Data Sheet

Product Name: Thalidomide-NH-CBP/p300 ligand 2

Cat. No.: GC63543

Chemical Properties

Cas. No. 2484739-21-9


Formula $C_{48}H_{57}F_2N_{11}O_6$ M.Wt 922.03

Solubility DMSO : 113.33 mg/mL (122.91 mM; Storage Store at -20°C, sealed storage, away
Need ultrasonic) from moisture and light

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

Thalidomide-NH-CBP/p300 ligand 2 (P-007) is a PROTAC-based CBP and p300 degrader (extracted from patent WO2020173440)[1].

[1]. LIU JING, et al. CYCLIC-AMP RESPONSE ELEMENT BINDING PROTEIN (CBP) AND/OR ADENOVIRAL E1A BINDING PROTEIN OF 300 KDA (P300) DEGRADATION COMPOUNDS AND METHODS OF USE. Patent WO2020173440.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA