

Product Data Sheet

Product Name: Begelomab

Cat. No.: GC74514

Chemical Properties

Cas. No. 1403744-56-8

Formula M.Wt

Solubility Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Begelomab

Background

Begelomab (SAND-26) is a murine IgG2b monoclonal antibody against DPP-4/CD26. Begelomab can be used for the research of severe refractory idiopathic inflammatory myopat.

References:

[1]. Lorenzo RD, et, al. Begelomab for severe refractory dermatomyositis: A case report. Medicine

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

(Baltimore). 2021 Mar 5;100(9):e24372.

[2]. Ojha R, et, al. Designing of a bispecific antibody against SARS-CoV-2 spike glycoprotein targeting human entry receptors DPP4 and ACE2. Hum Immunol. 2022 Apr;83(4):346-355.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA