

## Product Data Sheet

Product Name: H-D-Ala-D-Phe-OH hydrate

Cat. No.: GA22183

### Chemical Properties

Cas. No. 3061-94-7

Formula  $C_{12}H_{16}N_2O_3$

M.Wt 236.27

Solubility Soluble in DMSO

Storage Store at -20°C

**General tips** For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

**Shipping Condition** Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



### Background

Good substrate for D-peptidase S, an intracellular carboxypeptidase-like enzyme that preferentially hydrolyzes dipeptides containing hydrophobic D-amino acids.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA