

Product Data Sheet

Product Name: YAP-TEAD-IN-2

Cat. No.: GC73349

Chemical Properties

Cas. No. 2714432-83-2


Formula $C_{25}H_{24}ClFN_2O_4$ M.Wt 470.92

Solubility DMSO : 100 mg/mL (212.35 mM; Storage -20°C
Need ultrasonic)

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

YAP-TEAD-IN-2 (compound 6) is a potent YAP-TEAD PPI (protein-protein interaction) inhibitor with IC50 is 2.7 nM.

YAP-TEAD-IN-2

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

References:

[1]. Furet P, et al. The First Class of Small Molecules Potently Disrupting the YAP-TEAD Interaction by Direct Competition. ChemMedChem. 2022 Oct 6;17(19):e202200303.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA