

Product Data Sheet

Product Name: 2,2'-Biquinoline-4,4'-dicarboxylic acid dipotassium salt trihydrate

Cat. No.: GD03276

Chemical Properties

Cas. No. 207124-63-8

Formula C₂₀H₁₀K₂N₂O₄ · 3H₂O

M.Wt 474.55


Solubility

Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA