

Product Data Sheet

Product Name: Antiviral agent 55

Cat. No.: GC73805

Chemical Properties

Cas. No. 371137-60-9


Formula $C_{21}H_{20}N_2O_4$ M.Wt 364.39

Solubility DMSO : 3.7 mg/mL (10.15 mM; Need Storage -20°C
ultrasonic and warming)

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

Antiviral agent 55 (Compound 95) is an inhibitor for human immunodeficiency virus 1 and 2 (HIV 1 and HIV 2), and exhibits antiviral activity.

References:

[1]. Guenther RH, et al., Antiviral compounds and methods of use thereof. WO2011113060.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA