

Product Data Sheet

Product Name: 2H-Tetrazole,5-[4- (1H-1,2,4-triazol-1-yl) phenyl]

Cat. No.: GD14801

Chemical Properties

Cas. No. 1423810-54-1


Formula M.Wt

Solubility Storage argon filling

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA