

Product Data Sheet

Product Name: JF646-Hoechst

Cat. No.: GC71800

Chemical Properties

Cas. No. 2462011-60-3


Formula $C_{64}H_{70}N_{10}O_7Si$ M.Wt 1119.39

Solubility DMSO : 10 mg/mL (8.93 mM; Need ultrasonic and warming; Hygroscopic
Storage -20°C, protect from light, stored under nitrogen
DMSO has a significant impact on the solubility of product, please use newly opened DMSO)

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

JF646-Hoechst is a fluorescent red DNA probe that is an ideal substitute for large oligonucleotide-coupled antibodies used in PAINT experiments, especially for bacterial studies. JF646-Hoechst excitation/emission maximum =655/670 nm.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

References:

[1]. Spahn CK, et al. A toolbox for multiplexed super-resolution imaging of the E. coli nucleoid and membrane using novel PAINT labels. Sci Rep. 2018 Oct 3;8(1):14768.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA