

Product Data Sheet

Product Name: Thalidomide-NH-C9-NH₂ hydrochloride

Cat. No.: GC73617

Chemical Properties

Cas. No. 2305936-77-8

Formula C₂₂H₃₁ClN₄O₄


M.Wt 450.96

Solubility Storage -20°C, sealed storage, away from moisture and light

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Background

Thalidomide-NH-C9-NH₂ hydrochloride is the Thalidomide-based cereblon ligand used in the recruitment of CRBN protein. Thalidomide-NH-C9-NH₂ hydrochloride can be connected to the ligand for protein by a linker to form PROTACs.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA