

Product Data Sheet

Product Name: (S)-SKBG-1

Cat. No.: GC73510

Chemical Properties

Cas. No. 2955618-24-1


Formula $C_{22}H_{26}ClN_3O_6S$ M.Wt 495.98

Solubility DMSO : 100 mg/mL (201.62 mM; Storage -20°C
Need ultrasonic)

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

(S)-SKBG-1 is an inactive and covalent NONO ligand, used for assay control..

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA