

Product Data Sheet

Product Name: [Bicyclo[2.2.1]hept-5-en-2-yl]triethoxysilane (mixture of isomers)

Cat. No.: GD03775

Chemical Properties

Cas. No. 18401-43-9

Formula C₁₃H₂₄O₃Si

M.Wt 256.413


Solubility

Storage argon filling

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA