

Product Data Sheet

Product Name: Renditatug

Cat. No.: GC81436

Chemical Properties

Cas. No. 3037073-23-4

Formula M.Wt 145.42

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Renditatug

Background

Renditatug (BAT8008 Antibody) is an antibody targeting TROP2. Renditatug can be used for the synthesis of ADC BAT8008. The corresponding isotype control is: Human IgG1 kappa, Isotype Control [1] [2].

In Vitro, Renditatug is an ADC antibody targeting TROP2.

References:

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

[1]. Jianli Zhao et al. BAT8008, a TROP-2 antibody-drug conjugate (ADC), in patients with advanced solid tumor: Results from a phase 1 study.. J Clin Oncol 43, 3024-3024(2025).

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA