

Product Data Sheet

Product Name: DL-Tropic Acid

Cat. No.: GD21480

Chemical Properties

Cas. No. 552-63-6

Formula $C_9H_{10}O_3$

M.Wt 166.17


Solubility DMSO: 100mg/mL

Storage

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

Tropic acid (DL-Tropic acid) is a laboratory reagent used in the chemical synthesis of Atropine and Hyoscyamine.

[1]. Maria Gadzikowska, et al. Tropane Alkaloids in Pharmaceutical and Phytochemical Analysis. Acta Pol Pharm. Mar-Apr 2002;59(2):149-60.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA