

Product Data Sheet

Product Name: 2-Phenyl-4, 6-bis[4-(9, 9-dimethyl-9, 10-dihydroacridine)phenyl]-pyrimidine
Cat. No.: GF34299

Chemical Properties

Cas. No. 1870041-75-0

Formula C₅₂H₄₂N₄

M.Wt 722.94

Solubility Storage Store at RT, stored under nitrogen

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA