

Product Data Sheet

Product Name: PF-06795071

Cat. No.: GC70327

Chemical Properties

Cas. No. 2075629-81-9


Formula $C_{18}H_{17}F_4N_3O_3$ M.Wt 399.34

Solubility DMSO : 50 mg/mL (125.21 mM);
Need ultrasonic; Hygroscopic DMSO
has a significant impact on the
solubility of product, please use
newly opened DMSO)

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

PF-06795071 is a potent and selective covalent MAGL inhibitor with an IC₅₀ of 3 nM.

PF-06795071 shows excellent selectivity against FAAH with an IC₅₀ of 3.1 μM^[1].

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA