

Product Data Sheet

Product Name: Fesoterodine

Cat. No.: GC39379

Chemical Properties

Cas. No. 286930-02-7

SMILES O=C(C(C)C)OC1=CC=C(CO)C=C1[C@@H](C2=CC=CC=C2)CCN(C(C)C)C(C)C


Formula $C_{26}H_{37}NO_3$ M.Wt 411.58

Solubility Ethanol: 100 mg/mL (242.97 mM) Storage Store at -20°C, unstable in solution, ready to use.

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Background

Fesoterodine is an efficacious antimuscarinic agent for the overactive bladder (OAB)[1].

[1]. Ellsworth P, et al. Fesoterodine for the treatment of urinary incontinence and overactive bladder. Ther Clin Risk Manag. 2009;5:869-76. Epub 2009 Nov 18.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA