

Product Data Sheet

Product Name: Fmoc-Cys(Trt)-Wang resin

Cat. No.: GA10566

Chemical Properties

Cas. No. N/A

Chemical Name 9H-fluoren-9-ylmethyl N-[(2R)-1-oxo-3-tritylsulfanylpropan-2-yl]carbamate

SMILES C1=CC=C(C=C1)C(C2=CC=CC=C2)(C3=CC=CC=C3)SCC(C=O)NC(=O)OCC4C5=CC=CC=C5C6=CC=CC=C46


Formula 100-200:mesh,1%DVB M.Wt 0.3-0.8mmol/g

Solubility Storage Store at 2-8°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

Fmoc-Cys(Trt)-Wang resin

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA