

Product Data Sheet

Product Name: CCR2 antagonist 3

Cat. No.: GC31911

Chemical Properties

Cas. No. 1380100-86-6

SMILES CC1=C(F)C=CC(C(N(C)[C@@H](C(C)C)CN2CC(O)C2)=O)=C1


Formula $C_{17}H_{25}FN_2O_2$ M.Wt 308.39

Solubility Methanol : 250 mg/mL (810.66 mM; Need ultrasonic); DMSO : 100 mg/mL (324.26 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

CCR2 antagonist 3 is a chemokine receptor 2 (CCR2) antagonist.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA