

Product Data Sheet

Product Name: (R,S)-YAP-TEAD-IN-2

Cat. No.: GC81349

Chemical Properties

Cas. No.

Formula $C_{25}H_{24}ClFN_2O_4$

M.Wt 470.92


Solubility

Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

(R,S)-YAP-TEAD-IN-2 is the racemic mixture of YAP-TEAD-IN-2. YAP-TEAD-IN-2 (Compound 6) is a potent inhibitor of YAP-TEAD PPIs (protein-protein interaction) with an IC_{50} of 2.7 nM [1].

References:

[1]. Furet P, et al. The First Class of Small Molecules Potently Disrupting the YAP-TEAD Interaction by Direct Competition. ChemMedChem. 2022 Oct 6;17(19):e202200303.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA