

Product Data Sheet

Product Name: Boc-Gly-Gly-Phe-Gly-OH TFA

Cat. No.: GC35539

Chemical Properties

Cas. No. 2450273-39-7


Formula $C_{22}H_{29}F_3N_4O_9$ M.Wt 550.48

Solubility Water: 125 mg/mL (227.07 mM) Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

Boc-Gly-Gly-Phe-Gly-OH TFA is a self-assembly of N- and C-protected tetrapeptide.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA