

Product Data Sheet

Product Name: PROTAC EGFR degrader 7 diTFA

Cat. No.: GC73142

Chemical Properties

Cas. No.

Formula $C_{50}H_{50}F_6N_{10}O_{10}$

M.Wt 1064.98


Solubility DMSO : 100 mg/mL (93.90 mM;
Need ultrasonic)

Storage -20°C, sealed storage, away from
moisture and light

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

PROTAC EGFR degrader 7 (compound 13b) is a potent and selective CRBN-recruiting PROTAC EGFR L858R/T790M degrader, with a DC₅₀ of 13.2 nM. PROTAC EGFR degrader 7 inhibits NCI-H1975 cells proliferation, with an IC₅₀ of 46.82 nM. PROTAC EGFR degrader 7 significantly induces apoptosis and G2/M phase arrest in NCI-H1975 cell. PROTAC EGFR degrader 7 shows antitumor activity, and can be used for non-small cell lung cancer (NSCLC) research.

References:

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

[1]. Zhang W, et al. Discovery of highly potent and selective CRBN-recruiting EGFR^{L858R/T790M} degraders in vivo. Eur J Med Chem. 2022 Aug 5;238:114509.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA