

Product Data Sheet

Product Name: Cyclopropylboronic acid (contains varying amounts of Anhydride)

Cat. No.: GD06213

Chemical Properties

Cas. No. 411235-57-9

Formula

M.Wt

Solubility DMSO : 125 mg/mL (1455.18 mM; Storage -20°C argon filling
Need ultrasonic)

General For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the
tips ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Evaluation sample solution : ship with blue ice All other available size: ship with RT , or
Condition blue ice upon request.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA