

Product Data Sheet

Product Name: Copper(II) sulfate hydrate

Cat. No.: GB57864

Chemical Properties

Cas. No. 23254-43-5

Formula $\text{Cu.H}_2\text{O}_4\text{S.xH}_2\text{O}$

M.Wt 159.61 (anhydrous basis)

Solubility

Storage Store at room temperature, keep dry and cool

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA