

Product Data Sheet

Product Name: PROTAC Her3 Degradar-8

Cat. No.: GC73489

Chemical Properties

Cas. No. 2103331-95-7


Formula $C_{49}H_{55}N_{11}O_6S$ M.Wt 926.1

Solubility DMSO : 100 mg/mL (107.98 mM; Storage 4°C, protect from light
Need ultrasonic)

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

PROTAC Her3 Degradar-8 (Compound PP2) is a PROTAC that can degrade Her3 protein in vitro and cell experiments, and can be used to study diseases regulated by HER family proteins.

References:

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

[1]. Bifunctional molescules for her3 degradation and methods of use.WO2017117473.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA