

Product Data Sheet

Product Name: Phosphinous acid, bis(1, 1-dimethylethyl)-, 5-methoxy-1, 3-phenylene ester
Cat. No.: GF04705

Chemical Properties

Cas. No. 671235-32-8

Formula C₂₃H₄₂O₃P₂

M.Wt 428.53


Solubility

Storage Store at RT, sealed storage, away from moisture

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA