

Product Data Sheet

Product Name: Odronextamab

Cat. No.: GC72332

Chemical Properties

Cas. No. 1801338-64-6

Formula M.Wt

Solubility Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Odronextamab

Background

Odronextamab is a hinge-stabilised, fully human IgG4-based CD20 × CD3 bispecific antibody that binds CD3 on T cells and CD20 on B cells.

References:

[1]. Rajat Bannerji, et al. Odronextamab, a human CD20×CD3 bispecific antibody in patients with

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

CD20-positive B-cell malignancies (ELM-1): results from the relapsed or refractory non-Hodgkin lymphoma cohort in a single-arm, multicentre, phase 1 trial. Lancet H

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA