

Product Data Sheet

Product Name: SMARCA2/4-IN-2

Cat. No.: GC75800

Chemical Properties

Cas. No. 1339378-93-6


Formula $C_{10}H_9N_3O$ M.Wt 187.2

Solubility DMSO: 100 mg/mL (534.19 mM; Need ultrasonic) Storage $-20^{\circ}C$

General tips For obtaining a higher solubility , please warm the tube at $37^{\circ}C$ and shake it in the ultrasonic bath for a while. Stock solution can be stored below $-20^{\circ}C$ for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

SMARCA2/4-IN-2 is a PROTAC target protein ligand (Ligands for Target Protein for PROTACs).
SMARCA2/4-IN-2 can be used for synthesis PROTAC SMARCA2/4-degrader-5 [1].

References:

[1]. Preparation of heterocyclic compounds as SMARCA degraders and uses thereof. USA,
US20230072658A1 2020-12-17

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA