

Product Data Sheet

Product Name: STX-721

Cat. No.: GC71422

Chemical Properties

Cas. No.


Formula $C_{32}H_{35}ClN_6O_3$ M.Wt 587.11

Solubility DMSO : 100 mg/mL (170.33 mM; Storage -20°C
Need ultrasonic; Hygroscopic DMSO
has a significant impact on the
solubility of product, please use
newly opened DMSO)

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

STX-721 is a mutant-selective inhibitor of EGFR and HER2 Exon 20 insertion (ex20ins) mutants.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA