

Product Data Sheet

Product Name: 4-Methylanisole

Cat. No.: GC63288

Chemical Properties

Cas. No. 104-93-8


Formula $C_8H_{10}O$ M.Wt 122.16

Solubility DMSO : 100 mg/mL (818.60 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

4-Methylanisole (4-Methoxytoluene) is food flavoring agent and can be naturally found in Ylang Ylang fragrance oil[1].

[1]. Esther F A Brandon, et al. Implementation of Toxicokinetics in Toxicity studies--Toxicokinetics of 4-methylanisole and Its Metabolites in Juvenile and Adult Rats. Regul Toxicol Pharmacol. 2015 Oct;73(1):55-64.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA