

Product Data Sheet

Product Name: Fmoc-2,4-dimethoxy-4'-(carboxymethoxy)-benzhydramine

Cat. No.: GA21475

Chemical Properties

Cas. No. 126828-35-1

Formula $C_{32}H_{29}NO_7$

M.Wt 539.59

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

This Fmoc-amide linker (often called Knorr linker or Rink amide linker (RAM)) consists of a substituted benzhydramine residue attached to a moiety which can be coupled to commercially available amino-functionalized resins. It constitutes a very useful

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA