

Product Data Sheet

Product Name: Isopropylmagnesium chloride-lithium chloride complex solution(1.3 M in THF)
Cat. No.: GD15057

Chemical Properties

Cas. No. 745038-86-2


Formula $C_3H_7Cl_2LiMg$ M.Wt 145.24

Solubility Storage Store at RT, argon filling

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA