

Product Data Sheet

Product Name: KM04416

Cat. No.: GC69336

Chemical Properties

Cas. No. 26530-26-7

Formula $C_{12}H_{11}NO_3S$

M.Wt 249.29


Solubility DMSO : 50 mg/mL (200.57 mM; ultrasonic and warming and heat to 60°C)

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA