

Product Data Sheet

Product Name: MGR1
Cat. No.: GC45513

Chemical Properties

Cas. No. 2361529-46-4

SMILES O=C1C2C(C3C=CC2CC3)C(C4=C(OCOC(C(C)(C)C)=O)C=CC=C41)=O


Formula C₂₂H₂₄O₅ M.Wt 368.4

Solubility DMF: 30 mg/mL, DMSO: 30 mg/mL, Ethanol: 30 mg/mL
Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Background

MGR1 is a probe that generates reactive oxygen species (ROS).¹ Upon activation by esterases, MGR1 produces superoxide in vitro, an effect that is reversed by superoxide dismutase (SOD). MGR1 produces ROS and increases production of oxidized phosphatidylserine species in HEK293T cells. It reduces viability of HEK293T cells (IC₅₀s = 5.1-6.3 μM).

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

References

1. Kelkar, D., Ravikumar, G., Mehendale, N., et al. A chemical-genetic screen identifies ABHD12 as an oxidized-phosphatidylserine lipase. Nat. Chem. Biol. 15(2), 169-178 (2019).

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA