

Product Data Sheet

Product Name: (4-((3, 4-dihydroquinolin-1(2H)-yl)methyl)phenyl)boronic acid

Cat. No.: GF52149

Chemical Properties

Cas. No. 1029439-59-5

Formula C₁₆H₁₈BN₂O₂

M.Wt 267.13062

Solubility

Storage Store at RT

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA