

Product Data Sheet

Product Name: Gotistobart

Cat. No.: GC74431

Chemical Properties

Cas. No. 2226344-78-9

Formula

M.Wt

Solubility

Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Gotistobart

Background

Gotistobart (ONC-392) is a humanized anti-CTLA-4 antibody that confers immunotherapeutic effect by selective depletion of regulatory T cells (Treg) in the tumor microenvironment.

References:

[1]. Kai He, Safety and clinical activity of target-preserving anti-CTLA-4 antibody ONC-392 as

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

monotherapy in NSCLC patients who progressed on PD(L)1-targeted immunotherapy. Journal of Clinical Oncology, Volume 41, Issue 16_suppl.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA