

Product Data Sheet

Product Name: 2-methoxypyrazine

Cat. No.: GD16577

Chemical Properties

Cas. No. 3149-28-8


Formula M.Wt

Solubility DMSO : 100 mg/mL (908.18 mM; Storage
Need ultrasonic)

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

2-Methoxypyrazine is an active compound. 2-Methoxypyrazine can be used for the research of various biochemical studies.

[1]. C Sala, et al. Headspace solid-phase microextraction method for determining 3-alkyl-2-methoxypyrazines in musts by means of polydimethylsiloxane-divinylbenzene fibres. J Chromatogr A. 2000 Jun 2;880(1-2):93-9.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA