

Product Data Sheet

Product Name: SMCC Activatrd APC

Cat. No.: GC80378

Chemical Properties

Cas. No.

Formula

M.Wt

Solubility

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

SMCC Activatrd APC

Background

SMCC Activatrd APC is a pre-activated fluorescent labeling reagent (Ex/Em = 652/662 nm) that has been pre-modified with the heterodially crosslinking agent SMCC using alloxanthin protein

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

(Allophycocyanin, APC).

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA