

Product Data Sheet

Product Name: TRPC3 agonist-1

Cat. No.: GC79529

Chemical Properties

Cas. No. 1589567-51-0


Formula $C_{22}H_{23}F_3N_4O_3$ M.Wt 448.44

Solubility DMSO: 10 mg/mL (22.30 mM; Need warming) Storage Store at 4°C, stored under nitrogen

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Background

TRPC3 agonist-1 (compound 4n) is a selective and competitive TRPC3 channel agonist. TRPC3 agonist-1 binds at L2 site of hTRPC3 [1].

References:

[1]. Chen Y, et al. Structural mechanism of the agonist binding on human TRPC3 channel. Nat Commun. 2025 Oct 22;16(1):9343.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA