

Product Data Sheet

Product Name: Bis[1-Hydroxypyridine-2(1H)-Thionato-S,O]Copper(II)

Cat. No.: GD05297

Chemical Properties

Cas. No. 14915-37-8

Formula C₁₀H₈CuN₂O₂S₂

M.Wt 315.859


Solubility

Storage Store at room temperature away from moisture

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA