

Product Data Sheet

Product Name: Radanstobart

Cat. No.: GC79348

Chemical Properties

Cas. No. 3064483-12-8

Formula M.Wt 146.05

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Radanstobart

Background

Radanstobart (HFB-200603) is a humanized monoclonal antibody targeting CD272 with antitumor activity. Radanstobart is of the human IgG4k isotype, and the recommended isotype control is Human IgG4 (S228P) kappa, Isotype Control [1].

References:

[1]. WHO Drug Information, Vol. 39, No. 4, 2025

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA