

Product Data Sheet

Product Name: 2-Methyl-6-(piperidin-3-yl)pyridine-4-boronic acid

Cat. No.: GF34642

Chemical Properties

Cas. No. 2225153-02-4

Formula C₁₁H₁₇BN₂O₂

M.Wt 220.07588

Solubility

Storage -20°C, sealed storage, away from moisture

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA