

Product Data Sheet

Product Name: SERT/5-HT-IN-1

Cat. No.: GC79780

Chemical Properties

Cas. No. 3120534-54-2

Formula $C_{15}H_{18}N_2S$ M.Wt 258.38

Solubility DMSO: 25 mg/mL (96.76 mM; Need ultrasonic) Storage Store at 4°C, away from moisture

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.


Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Structure



Background

SERT/5-HT-IN-1 (compound D17) is an orally active and potent dual target inhibitor of SERT/5-HT3A (SERT IC₅₀ = 1.5 nM, NET IC₅₀ = 49 nM, DAT IC₅₀ = 91 nM, 5-HT3A K_i = 12 nM). SERT/5-HT-IN-1 has an IC₅₀ of 39 nM to 5-HT3A. SERT/5-HT-IN-1 has a lower risk of liver, kidney, and cardiac toxicity. SERT/5-HT-IN-1 can be used for research on depressive disorders [1].

References:

[1]. Zheng J, et al. Design, synthesis, and biological profiling of novel aryl-spirocyclic diamine derivatives with potential antidepressant-like properties. Eur J Med Chem. 2026 Jan 5;301:118198.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA