

## Product Data Sheet

Product Name: LITHIUM 2-THIENYLCYANOCUPRATE

Cat. No.: GF10297

### Chemical Properties

Cas. No. 112426-02-5

Formula  $C_5H_3CuLiNS$

M.Wt 179.64

Solubility

Storage 2-8°C, sealed storage, away from moisture,

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA