

Product Data Sheet

Product Name: PS432
Cat. No.: GC72880

Chemical Properties

Cas. No. 2083630-26-4


Formula $C_{25}H_{19}ClN_2O_5S$ M.Wt 494.95

Solubility DMSO : 100 mg/mL (202.04 mM; Need ultrasonic) Storage $-20^{\circ}C$

General tips For obtaining a higher solubility , please warm the tube at $37^{\circ}C$ and shake it in the ultrasonic bath for a while. Stock solution can be stored below $-20^{\circ}C$ for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

PS432 is a PKC inhibitor with IC₅₀s of 16.9 μM (PKC α) and 18.5 μM (PKC ζ), respectively. PS432 effectively inhibits the proliferation of non-small cell lung cancer cells (NSCLCs) and tumor growth in mouse xenograft models.

PS432 (25 μM ; 24 h) inhibits cancer cells proliferation, with IC₅₀s of 30.7 μM (A549) and 12.4 μM

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

(A427) in semi solid agar medium, for example[1]. PS432 (50 μ M; 12, 24, and 36 h) arrests cell cycle at G0/G1 phase in A549 lung cancer cells[1].

PS432 (2.5 mg/kg; i.p.; once daily for 14 d) inhibits tumor growth in xenograft model of lung cancer[1].

References:

[1]. Arencibia JM, et al. An Allosteric Inhibitor Scaffold Targeting the PIF-Pocket of Atypical Protein Kinase C Isoforms. ACS Chem Biol. 2017 Feb 17;12(2):564-573.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA