

Product Data Sheet

Product Name: (S)-(+)-3-Fluoropyrrolidine hydrochloride

Cat. No.: GD15430

Chemical Properties

Cas. No. 136725-53-6

Formula

M.Wt


Solubility

Storage 2-8°C argon filling

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

(S)-3-Fluoropyrrolidine hydrochloride is a biochemical reagent that can be used as a biological material or organic compound for life science related research.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA